

4.5" x 9.5" Envelope

Size 4.5" x 9.5"  
(114,3mm x 241.3mm)

Suggested Paper Stock:  
Bagas Paper 90kg  
\* if not available, recycled paper

Colors for Printing

 KYOCERA Red  
( Standardized color: C5% + M100% + Y90% + K0%  
Closest color: PANTONE 1797C )

 KYOCERA Black ( Black 100% )

\* The official corporate logotype  
can be downloaded from the  
"Download" section of this  
website.

\* The width of the Brand Symbol  
should be 0.2A, where A equals  
the length of the envelope.

Shown actual size

Unit : mm

